



<b>Description:</b>	Rescue sheet
<b>Final dimensions:</b>	2.1 m x150 cm
<b>Composition:</b>	Polypropylene.
<b>Latex:</b>	Free
<b>PVC:</b>	Free
<b>Re use:</b>	No
<b>Sterilizable by:</b>	Ethylene oxide
<b>Thickness:</b>	0,015±10% mm
<b>Colour:</b>	Silver-Gold
<b>Shelf life:</b>	120 months
<b>CFU/g</b>	<10 <sup>2</sup>



<b>Nr and Revision Date:</b>	8 -- 2019/04/30
<b>Directive:</b>	Not applicable
<b>Applicable norms:</b>	EN ISO 13485:2016 ; EN ISO 11737-1:2006 ; ENISO11737-1:2006/AC:2009 ;

Level	Type	Dimension	Stamp	Labels
1	1 Flexible peel-open pack	134x200mm		labeled
2	100 <b>Carton box</b>	<b>769x381x110mm</b>		<b>labeled</b>

<b>Storage instructions</b>	Store in a dry place at room temperature and protect from sunlight. The product has good stability and therefore, normal conditions of medical devices storage shall apply.
<b>Waste instructions</b>	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
<b>Intended use</b>	For protection in emergencies

<b>WASTE AND ENVIRONMENT REGULATIONS</b>	Weight (gr) of packing material used in 100 pcs	1452
--	---	------

**ISSUED BY: Quality Management Department, 2019/08/08**

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.