

COLPOSCOPE OP-C2

THE BEST VISION IN COLPOSCOPY

- 3 different magnifications
- Easy to handle
- Excellent illumination



OPTOMIC[®]
producing quality

The OPTOMIC OP-C2 colposcope is the perfect equipment for diagnostic observations as well as minor routine surgeries in the doctor's practice, helping the specialist in his everyday work.

COLPOSCOPE **OP-C2**



INCLINED BINOCULAR HEAD



REMOVABLE COLD LIGHT SOURCE



OPTIONAL VIDEO SYSTEM

- 45° inclined stereoscopic binocular head.
- One central handle to operate the head unit easily.
- The arm is an autocompensated parallelogram which can be placed in any position.
- Very soft and smooth operation of all functions.
- Cold light sources of Halogen, Xenon or LED light; all with swing-in green filter.
- Removable cold light source appropriate for other applications such as endoscopy, illumination for speculum or head lamp.
- High intensity fiber optic illumination.
- Wide field eyepieces with dioptic adjustment and protective rubber shells.
- 3 steps Galilei magnifications changer and micrometric focusing system.
- Working distance of 250mm with standard objective (f=300mm optional).
- Range of magnifications: 3.5x up to 17.5x, depending on the chosen combination of objective, binocular and eyepieces.
- Different floor and wall mounts available.
- Optional monitor supports and trays for laptop or CCD camera.
- Optional systems for video-colposcopy or digital photography.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
export@optomic.com • www.optomic.com



producing quality

Specifications are subject to change without notice

