



Veterinary scale "Cosy"

8310.02

- ▶ Memory function to recall the last weight value
- ▶ Easy to clean plastic housing
- ▶ Hold and tare function
- ▶ Energy-saving auto-off function

Features

- ▶ Simple operation
- ▶ Easy to read LC display, digit height 25 mm
- ▶ Easy to clean plastic housing
- ▶ Tare function: Covers/paddings can be tared
- ▶ Hold function: The weight value is frozen, after unloading, the scale auto-matically returns to zero
- ▶ Memory function: The last weight value determined can be recalled
- ▶ Units in kg or lb
- ▶ Energy-saving auto-off function
- ▶ Battery operation with four AA 1.5 V batteries (included)
- ▶ Power supply operation optional (see accessories)

Products

Order number Variation	Approvable	Division g	Maximum load kg
8310.02.001 Baby scale "Cosy"	No	5	20

Accessories

Order number
Variation

Description

5908.09.002

Power supply 6V / 500 mA
