


vyaire™
MEDICAL




AirLife™ Misty Fast™
small-volume nebulizer

Support your patient's next breath while reducing treatment time



AirLife™ Misty Fast™ disposable SVN

The AirLife Misty Fast™ small-volume nebulizer (SVN) is a small volume, high-efficiency air entrainment nebulizer that delivers comparable results to the AirLife Misty Max 10™ nebulizer, but in half the time. The one-piece jet design provides efficient, high output rates for fast treatment times and consistent medication delivery. Nine kit configurations accommodate your patient's aerosolization needs, and all kits include our crush-resistant O₂ tubing.



The AirLife Misty Fast SVN offers the following key benefits:

- Treatment requires half the time of the AirLife Misty Max 10 SVN
- Proprietary swivel mouthpiece helps enhance patient positioning
- Unique, integrated jet and baffle design provides medication delivery consistency and reliability
- Clear base with visible graduations facilitates medication loading and aerosol delivery visualization

To view a product demonstration or calculate how much time that you could save with MistyFast, please visit us at vyaire.com/MistyFast.

AirLife Misty Fast SVN (50 per case)

Cat. no.	Product description
FN2431	Misty Fast SVN and 7' O ₂ tubing
FN2433	Misty Fast SVN with adult aerosol mask and 7' O ₂ tubing
FN2434	Misty Fast SVN with baffled tee swivel mouthpiece and 7' O ₂ tubing
FN2438	Misty Fast SVN with baffled tee swivel mouthpiece, 6" flextube and 7' O ₂ tubing
FN2439	Misty Fast SVN with baffled tee swivel mouthpiece, 6" flextube and 14' U/Connect-It® O ₂ tubing
FN2444	Misty Fast SVN with pediatric aerosol mask and 7' O ₂ tubing
FN2445	Misty Fast SVN with baffled tee swivel mouthpiece and 7' U/Connect-It O ₂ tubing
FN2446	Misty Fast SVN with baffled tee swivel mouthpiece, 6" flextube and 7' U/Connect-It O ₂ tubing
FN2455	Misty Fast SVN with pediatric aerosol elephant mask and 7' O ₂ tubing



For more information, please contact your local Vyairē representative at **833.327.3284**.

GLOBAL
HEADQUARTERS

Vyairē Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045, USA

vyaire.com

Rx ONLY For U.S. distribution only.

© 2018 Vyairē. Vyairē, the Vyairē logo and all other trademarks or registered trademarks are property of Vyairē Medical, Inc., or one of its affiliates. VYR-US-1800012 (1118)