

# Surgical Suction

# Su


Suction Connecting Tubing  
Link Suction Connecting Tubing Suction Devices  
Poole Suckers  
Suction Probes  
Insufflation  
Single-use Medical Devices  
Uterine Aspiration Curette Systems  
Suction Collection Systems  
Accessories

# 02.

## Contents

Page

<b>Suction Connecting Tubing</b>	Suction Connecting Tubing	3
	Soft Suction Connecting Tubing	4
	Conductive Tubing with Anti-static Line	4
	Flexi-Rib Tubing	4
	Tubing with Capped Vacuum Control	5
	Capped Vacuum Control Connector	5
	5-in-1 Connector	5
<b>Link Suction Connecting Tubing</b>	Link Suction Connecting Tubing	6
	Box Link Suction Connecting Tubing	7
<b>Suction Devices</b>	Link Yankauers	8
	Link Yankauer Sets	9
	Derby Yankauers	10
	Derby Yankauer Sets	10
	Selectovac Yankauer	10
	Poole Suckers	11
<b>Poole Suckers</b>	Poole Sucker Sets	11
	Poole Suction & Flexi-Rib Tubing Sets	11
	Suction Probes	12
	Suction Probe Sets	12
<b>Suction Probes</b>	Micro-Flex and Neuro-Flex Suction Probes	13
	Endo-Flex Suction Devices	13
	Sigmo-Flex Suction Devices	13
	Frazier Probes	13
<b>Insufflation</b>	Insufflation Tubing	14
<b>Single-use Medical Devices</b>	Easyflow, Haemo-Ligator & Atiomo DyeSeal	15
<b>Uterine Aspiration Curette System</b>	Aspiration Curettes	16
	Karman Curettes	16
	Wide Bore Aspiration Tubing (Flexi-Cuff)	17
<b>Hy-G-Vac Suction Collection System</b>	Hy-G-Vac Suction Collection System	17
<b>Accessories</b>	Sharps Pad	18
	Tip Cleaner	18
	Diathermy Pencil	18
	Smoke Evacuation Diathermy Pencil Set	18
	Camera Drape	19
	C-Arm Drape	19



## Suction Connecting Tubing

A complete system to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices. Pennine offer a range of easily connected, interchangeable products in a choice of tubing length and diameter.


Sterile double wrapped individual peel pouches.

'F' (Female) connector fits suction machines and male tapered connections for instrumentation.

'M' (Male) connector push fits into the 'F' connector and is normally used to connect to the standard funnel fitted on suction catheters. Can also be used to connect lengths of tubing together for longer lengths.

Application of 'M' connector with 'VC' (vacuum control) vent converts a standard catheter to vacuum control application.


- F Female
- FF Female/Female Connectors
- FFM Female/Female/Male Connectors
- FFM/VC Female/Female/Male Connectors with Vacuum Control



Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
CT-4312	F/F	1.8m	24CH	5	25
CT-4362	F/F	2.7m	24CH	5	20
CT-4372	F/F	3.5m	24CH	5	20
CT-4314	F/F/M	1.8m	24CH	5	25
CT-4364	F/F/M	2.7m	24CH	5	20
CT-4334	F/F/M	3.0m	24CH	5	20
CT-4394/VC	F/F/M/VC	2.7m	24CH	5	20
CT-4374/VC	F/F/M/VC	3.5m	24CH	5	20
CT-4314/VC	F/F/M/VC	1.8m	24CH	5	25
CT-4324/VC	F/F/M/VC	2.0m	24CH	5	25
CT-4334/VC	F/F/M/VC	3.0m	24CH	5	20
CT-4012	F/F	1.8m	26CH	6	25
CT-4022	F/F	2.0m	26CH	6	25
CT-4042	F/F	2.7m	26CH	6	20
CT-4032	F/F	3.0m	26CH	6	20
CT-4072	F/F	3.5m	26CH	6	20
CT-4092	F/F	4.0m	26CH	6	15
CT-4062	F/F	6.0m	26CH	6	10
CT-4020	F/open end	2.0m	26CH	6	25
CT-4030	F/open end	3.0m	26CH	6	20
CT-4014	F/F/M	1.8m	26CH	6	25
CT-4024	F/F/M	2.0m	26CH	6	25
CT-4044	F/F/M	2.7m	26CH	6	20
CT-4034	F/F/M	3.0m	26CH	6	20
CT-4074	F/F/M	3.5m	26CH	6	20
CT-4064	F/F/M	6.0m	26CH	6	10
CT-4014/VC	F/F/M/VC	1.8m	26CH	6	25
CT-4024/VC	F/F/M/VC	2.0m	26CH	6	25
CT-4034/VC	F/F/M/VC	3.0m	26CH	6	20
CT-4622	F/F	2.0m	30CH	8	10
CT-4632	F/F	3.0m	30CH	8	10
CT-4672	F/F	3.5m	30CH	8	10
CT-4642	F/F	4.0m	30CH	8	10
CT-4624	F/F/M	2.0m	30CH	8	10
CT-4634	F/F/M	3.0m	30CH	8	10
CT-4674	F/F/M	3.5m	30CH	8	10
CT-4634/VC	F/F/M/VC	3.0m	30CH	8	10

Please note: All dimensions are nominal.

# 04.



- Softer PVC grade tubing for enhanced durability
- Not to be used in high pressure applications
- Maximum pressure 300mm

## Soft Suction Connecting Tubing


Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
CTS-4332	F/F	3.0m	24CH	5	20
CTS-4022	F/F	2.0m	26CH	6	25
CTS-4032	F/F	3.0m	26CH	6	20
CTS-4072	F/F	3.5m	26CH	6	20
CTS-4092	F/F	4.0m	26CH	6	15
CTS-4640/R	F/open end	4.0m	30CH	8	10



- Anti-static line, prevents discharge
- Standard PVC grade tubing

## Conductive Tubing with Anti-static Line


Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
CCT-4112	F/F	1.8m	26CH	6	25
CCT-4122	F/F	2.0m	26CH	6	25
CCT-4142	F/F	2.7m	26CH	6	25
CCT-4132	F/F	3.0m	26CH	6	20
CCT-4172	F/F	3.5m	26CH	6	20
CCT-4134	F/F/M	3.0m	24CH	6	20
CCT-4334	F/F/M	3.0m	26CH	5	20
CCT-4734	F/F/M	3.0m	28CH	7	50




- Ribbed connecting tubing for enhanced drapability
- Non-slip outer rib reduces friction
- Resistant to coiling, kinking
- Maximum pressure 500mm HG Mercury
- Colour coded funnel identification
- Multi-fit to all suction devices and sources
- Economical box size

## Flexi-Rib Tubing

Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
CT-4922	F/F	2.0m	26CH	6	60
CT-4932	F/F	3.0m	26CH	6	50
CT-4972	F/F	3.5m	26CH	6	50
CT-4924	F/F/M	2.0m	26CH	6	60
CT-4934	F/F/M	3.0m	26CH	6	50
CT-4722	F/F	2.0m	28CH	7	60
CT-4732	F/F	3.0m	28CH	7	50
CT-4772	F/F	3.5m	28CH	7	40
CT-4724	F/F/M	2.0m	28CH	7	60
CT-4734	F/F/M	3.0m	28CH	7	50
CT-4774	F/F/M	3.5m	28CH	7	40
CT-4714/VC	F/F/M/VC	1.8m	28CH	7	25
CT-4741/VC	F/F/M/VC	2.5m	28CH	7	20





## Tubing with Capped Vacuum Control

Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
CT-4015	F/VCU2	1.8m	26CH	6	25
CT-4025	F/VCU2	2.0m	26CH	6	25



## Capped Vacuum Control Connector

Fingertip vacuum control connection capped for ease of use to convert a standard funnel suction catheter to a vacuum controlled suction catheter.

- Available sterile and bulk non-sterile

Product Code	Type	Qty. Per Box
VCU-2	Sterile	100
VC-2/Bulk	Bulk, non Sterile	1200




## 5-in-1 Connector

- Clear PVC, single use only
- Can be used to extend length of tubing i.e. to connect two tubes together of either the same or different diameters
- Available sterile and bulk non-sterile
- Designed for use with all types of tubing

Product Code	Type	Qty. Per Box
DP-501	Sterile	100
DP-501/Bulk	Bulk, non Sterile	500

# 06.



## Link Suction Connecting Tubing / Bubble Tubing

The Pennine Link system is used for connecting suction sources to suction waste collection systems, suction catheters, yankauer suckers, suction probes and other suction devices. Tubing extruded with an integral female funnel and a tapered male end to provide adaptable connections. The tube can be modified by cutting with scissors.

This unique system will guarantee an easy and faultless connection between all the individual products in the Link range. Link Tubing is also referred to as Bubble Tubing.

Although the Link products are optimised to be used within the Link system, they can also be easily connected to other products in our range.

- Sterile, double wrapped transparent individual peel pouches
- Disposable PVC tubing
- Available in two diameters, 5mm & 7mm

Available either in plain PVC or with an antistatic black line to prevent static discharge (also known as conductive tubing)


Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
CT-4513	Plain F/F	1.8m	24CH	5mm	25
CT-4543	Plain F/F	2.7m	24CH	5mm	20
CT-4511	Plain M/F	1.8m	24CH	5mm	25
CT-4541	Plain M/F	2.7m	24CH	5mm	20
CT-4713	Plain F/F	1.8m	28CH	7mm	25
CT-4783	Plain F/F	3.7m	28CH	7mm	20
CT-4741	Plain M/F	2.7m	28CH	7mm	20
CT-4743	Plain F/F	2.7m	28CH	7mm	20
CT-4711	Plain M/F	1.8m	28CH	7mm	25
CCT-4513	Conductive F/F	1.8m	24CH	5mm	25
CCT-4543	Conductive F/F	2.7m	24CH	5mm	20
CCT-4583	Conductive F/F	3.7m	24CH	5mm	15
CCT-4512	Conductive M/F	1.8m	24CH	5mm	25
CCT-4542	Conductive M/F	2.7m	24CH	5mm	20
CCT-4582	Conductive M/F	3.7m	24CH	5mm	15
CCT-4713	Conductive F/F	1.8m	28CH	7mm	25
CCT-4743	Conductive F/F	2.7m	28CH	7mm	20
CCT-4783	Conductive F/F	3.7m	28CH	7mm	10
CCT-4712	Conductive M/F	1.8m	28CH	7mm	25
CCT-4742	Conductive M/F	2.7m	28CH	7mm	20
CCT-4782	Conductive M/F	3.7m	28CH	7mm	10



## Box Link Suction Connecting Tubing / Bubble Tubing

LINK tubing available in 30 and 50 metre lengths, three diameters: 5mm, 6mm and 7mm supplied in convenient dispenser boxes. Tube extruded with a bubble at every 0.90m intervals. Simply cut through the bubble at the required length to create either a male or female end which will fit a suction tip or suction catheter.

Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
SL-8005	Plain	30m	24CH	5mm	1 Roll
SL-8055	Plain	50m	24CH	5mm	1 Roll
SL-8056	Plain	50m	26CH	6mm	1 Roll
SL-8007	Plain	30m	28CH	7mm	1 Roll
SL-8205	Conductive	30m	24CH	5mm	1 Roll
SL-8255	Conductive	50m	24CH	5mm	1 Roll
SL-8256	Conductive	50m	26CH	6mm	1 Roll
SL-8207	Conductive	30m	28CH	7mm	1 Roll



For integral funnel connector - cut here



For integral funnel connector - cut here

# 08.



## Link Yankauers

Rigid, transparent suction devices for use with suction connecting tubing. For general purpose surgical suctioning. Also available as a set, pre-attached to 2 or 3 metre tubing.

- Sterile, double wrapped transparent individual peel pouches
- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip
- Midi Yankauer available with additional Atrauma Tip to reduce tissue grab



Short Mini (no side eyes)    Short Mini Super Tip (no side eyes)    Mini Fine Tip (two eyes)



Mini Fine Super Tip (two eyes)    Midi (six eyes)    Midi Super Tip (six eyes)



Maxi (six eyes)    Maxi Super Tip (six eyes)    Maxi Plus (four eyes)



Long Rigid




Short Rigid



Soft Rigid

Product Code	Type	Use	Length	CH	Qty. Per Box
LY-3400	Short Mini	Paediatric/dental	190mm	15ch	100
LY-3401	Short Mini VC		190mm	15ch	100
LY-4400	Short Mini Supertip		195mm	15ch	100
LY-4401	Short Mini Supertip VC		195mm	15ch	100
LY-3500	Mini Fine Tip	Fine suctioning (e.g. ENT)	270mm	12ch	50
LY-3501	Mini Fine Tip VC		270mm	12ch	50
LY-4500	Mini Fine Supertip		275mm	12ch	50
LY-4501	Mini Fine Supertip VC		275mm	12ch	50
LY-3600	Midi 30° angle	For general and pharyngeal suction	280mm	18ch	50
LY-3601	Midi VC		280mm	18ch	50
LY-3600/PR	Midi 10° angle		270mm	18ch	50
LY-4600	Midi Supertip		280mm	18ch	50
LY-4601	Midi Supertip VC		280mm	18ch	50
LY-3700	Maxi		For large particles and significant volumes of fluid	275mm	22ch
LY-3701	Maxi VC	275mm		22ch	50
LY-4700	Maxi Supertip	275mm		22ch	50
LY-4701	Maxi Supertip VC	275mm		22ch	50
LY-3800	Maxi Plus	For orthopaedic use	280mm	28ch	50
LY-3801	Maxi Plus VC		280mm	28ch	50
LY-3602	Atrauma tip, short soft	For general suction with minimal tissue grab	280mm	18ch	50
LY-3603	As above with VC		280mm	18ch	50
LY-3606	Atrauma tip, short rigid	For general suction with minimal tissue grab	280mm	18ch	50
LY-3607	As above with VC		280mm	18ch	50
LY-3608	Atrauma tip, long rigid	For general suction with minimal tissue grab	280mm	18ch	50
LY-3609	As above with VC		280mm	18ch	50





- Sterile, double wrapped transparent individual peel pouches
- Available with or without VC (vacuum control) vent

### Link Yankauers Sets

Link yankauers pre-connected to Pennine suction connection tubing in 2.0 and 3.0 metre lengths. For general purpose surgical suctioning.


#### Sets with 2m tubing

Code No.	Tubing Diam./Length	Qty. Per Box	Yankauer Type
LYS-5420	6mm / 2m	10	LY-3400 Short Mini
LYS-5421	6mm / 2m	10	LY-3401 Short Mini VC
LYS-5520	6mm / 2m	10	LY-3500 Mini Fine Tip
LYS-5521	6mm / 2m	10	LY-3500 Mini Fine Tip VC
LYS-5620	6mm / 2m	10	LY-3600 Midi 30° angle
LYS-5621	6mm / 2m	10	LY-3601 Midi VC
LYS-5720	6mm / 2m	10	LY-3700 Maxi
LYS-5721	6mm / 2m	10	LY-3701 Maxi VC
LYS-6820	6mm / 2m	10	LY-3800 Maxi Plus
LYS-6821	6mm / 2m	10	LY-3801 Maxi Plus VC
LYS-3622	6mm / 2m	10	LY-3602 Atrauma tip, short soft
LYS-3623	6mm / 2m	10	LY-3603 Atrauma tip, short soft VC
LYS-3626	6mm / 2m	10	LY-3606 Atrauma tip, short rigid
LYS-3627	6mm / 2m	10	LY-3607 Atrauma tip, short rigid VC
LYS-3628	6mm / 2m	10	LY-3608 Atrauma tip, long rigid
LYS-3629	6mm / 2m	10	LY-3609 Atrauma tip, long rigid VC

#### Sets with 3m tubing

Code No.	Tubing Diam./Length	Qty. Per Box	Yankauer Type
LYS-5430	6mm / 3m	10	LY-3400 Short Mini
LYS-5431	6mm / 3m	10	LY-3401 Short Mini VC
LYS-5530	6mm / 3m	10	LY-3500 Mini Fine Tip
LYS-5531	6mm / 3m	10	LY-3501 Mini Fine Tip VC
LYS-5630	6mm / 3m	10	LY-3600 Midi 30° angle
LYS-5631	6mm / 3m	10	LY-3601 Midi VC
LYS-5730	6mm / 3m	10	LY-3700 Maxi
LYS-5731	6mm / 3m	10	LY-3701 Maxi VC
LYS-6830	6mm / 3m	10	LY-3800 Maxi Plus
LYS-6831	6mm / 3m	10	LY-3801 Maxi Plus VC
LYS-3632	6mm / 3m	10	LY-3602 Atrauma tip, short soft
LYS-3633	6mm / 3m	10	LY-3603 Atrauma tip, short soft VC
LYS-3636	6mm / 3m	10	LY-3606 Atrauma tip, short rigid
LYS-3637	6mm / 3m	10	LY-3607 Atrauma tip, short rigid VC
LYS-3638	6mm / 3m	10	LY-3608 Atrauma tip, long rigid
LYS-3639	6mm / 3m	10	LY-3609 Atrauma tip, long rigid VC

# 10.



- Sterile, double wrapped in individual peel pouches
- Available with or without VC (Vacuum control) vent
- Available with or without sump tip to reduce the risk of tissue damage

## Derby Yankauers

Clear, acrylic suction devices, similar to traditional yankauer suckers. With or without sump tip to reduce potential tissue suction during use.

Product Code	Type	Length	CH	Inside Diam. mm	Qty. Per Box
YS-3004	with sump tip	265mm	24CH	8	50
YS-3005	with sump tip and VC	265mm	24CH	8	50
YS-3014	straight tip	265mm	24CH	8	50
YS-3015	straight tip and VC	265mm	24CH	8	50
YS-3000	Midi Yankauer with straight tip	265mm	20CH	6.66	50
YS-3001	Midi Yankauer with VC	265mm	20CH	6.66	50

## Derby Yankauers Sets

Derby Yankauer Suckers pre-attached to 2m or 3m connecting tube. Available with or without 'sump' tip to reduce the risk of tissue damage during use.

- All Pennine suction sets come sterile, double wrapped in individual peel pouches
- Available with or without VC (vacuum control) hole

Set with 2m tubing	CH/Tip	Qty. Per Box	Set with 3m tubing	Qty. Per Box	Tube Diam. mm	Yankauer type
SYS-5820	24CH	20	SYS-5830	10	6	YS-3004
SYS-5821	24CH	20	SYS-5831	10	6	YS-3005
SYS-5520	24CH	20	SYS-5530	15	6	YS-3014
SYS-5521	24CH	20	SYS-5531	15	6	YS-3015



- Switch system allows the operator to manually open and close suction
- Easy connection to Pennine range of tubing

## Selectavac Yankauer

Surgical suction cannula with switch on and off valve, available with bulb tip and straight tip.

Product Code	Type	Length	CH/Tip	Inside Diam. mm	Qty. Per Box
YS-3009	Bulb tip	265	24CH	8	50
YS-3010	Straight tip	265	24CH	8	50



Flexible straight tip



Flexible angled tip



Rigid tip

## Pooler Suckers


Flexible or Rigid sump sucker with removable inner suction probe. Specifically for abdominal surgery, reducing trauma to surrounding tissues.

- Available with atraumatic end
- Available with angled tip
- Removable inner probe

Product Code	Type	Length	CH/Tip	Inside Diam. mm	Qty. Per Box
YS-3006	Flexible straight tip	300	20CH	6.66	50
YS-3008	Flexible (angled tip)	300	20CH	6.66	50
YS-3007	Rigid (White)	300	15CH	6.5	50

## Pooler Sucker Sets

Flexible and rigid 'sump' suckers pre-attached to 2m or 3m connecting tube. For deep abdominal suctioning or where a 'sump' suction action is required.



Set with 2m tubing	CH/Tip	Qty. Per Box	Set with 3m tubing	Qty. Per Box	Tube Diam. mm	Pooler type
SYS-5920	20CH	10	SYS-5930	10	6	YS-3006
SYS-5923	20CH	15	SYS-5933	10	6	YS-3008
SYS-5720	15CH	10	SYS-5730	10	6	YS-3007

## Pooler Suction & Flexi-Rib Tubing Sets

3m Flexi-Rib tubing pre-attached to Pooler suction probe.

Set with 3m tubing	CH/Tip	Qty. Per Box	Tube Diam. mm	Pooler type
SYS-5931	20CH	50	6mm Flexi-Rib	YS-0006
SYS-5932	20CH	50	6mm Flexi-Rib	YS-0008

# 12.

SP-3800 - 4 side eyes



MP-3800 - 4 side eyes



SP-3900 - No side eyes




SP-4200 - No side eyes



SP-4300 - No side eyes



- Standard probes 20ch
- Available in a wider bore, 30ch for orthopaedic use and for use with wide bore suction connection tubing



## Suction Probes

Pennine suction devices manufactured from rigid transparent PVC, for general suction purposes, orthopaedic suctioning or dental work. Easily attached to Pennine suction tubing.


Product Code	Type	Length	CH/Tip	Qty. Per Box
SP-3800	Suction probe with male connector	245mm	20CH	50
SP-3801	As above with VC	245mm	20CH	50
SP-3900	Flat tip probe with male connector	250mm	20CH	50
SP-3901	As above with VC	250mm	20CH	50
SP-4200	Maxi probe	280mm	30CH	50
SP-4201	As above with VC	280mm	30CH	50
SP-4300	Flat tip Maxi probe	255mm	30CH	50
SP-4301	As above with VC	255mm	30CH	50
MP-3800	Dental suction device	120mm	20CH	50

## Suction Probe Sets

Suction probes available as sets pre-attached to 2m or 3m tubing. For general use (by attachment to 6mm diameter standard connecting tubing) or orthopaedic use (by attachment to 30ch, 8mm wide bore tubing).

- Standard probes 20ch
- Available in a wider bore, 30ch for orthopaedic use and for use with wide bore suction connection tubing

Set with 2m tubing	CH/Tip	Qty. Per Box	Tube Diameter	Probe type	Set with 2m tubing	CH/Tip	Qty. Per Box	Tube Diameter	Probe type
SPS-5020	20CH	25	6mm	SP-3800	SPS-5030	20CH	20	6mm	SP-3800
SPS-5021	20CH	25	6mm	SP-3801	SPS-5031	20CH	20	6mm	SP-3801
SPS-5320/G	20CH	25	6mm	SP-3900	SPS-5330/G	20CH	20	6mm	SP-3900
SPS-5321/G	20CH	25	6mm	SP-3901	SPS-5331	20CH	20	6mm	SP-3901
SPS-6020	30CH	10	8mm	SP-4200	SPS-6030	30CH	10	8mm	SP-4200
SPS-6021	30CH	10	8mm	SP-4201	SPS-6031	30CH	10	8mm	SP-4201
SPS-6320	30CH	10	8mm	SP-4300	SPS-6330	30CH	10	8mm	SP-4300
SPS-6321	30CH	10	8mm	SP-4301	SPS-6331	30CH	10	8mm	SP-4301




- Sterile, double wrapped individual peel pouches.

## Micro-Flex and Neuro-Flex Suction Probes

Special PVC disposable probes with efficient prestrol vacuum control device. Designed for neuro surgery, ENT surgery and in applications where controlled suction is required.

Product Code	Use	Length	CH/ Tip	Qty. Per Box
YS-3148/14	For neuro surgery and other intricate suction procedures.	16cm	14ch ø 4.66mm	50
YS-3128		15cm	12ch ø 4.00mm	50
YS-3108	For neuro surgery and other intricate suction procedures.	15cm	10ch ø 3.33mm	50
YS-3088		15cm	8ch ø 2.66mm	50
YS-3088(Z)	For otological procedures and other intricate surgical applications	11cm	8ch ø 2.66mm	50
YS-3068		11cm	6ch ø 2.00mm	50
YS-3058		11cm	5ch ø 1.66mm	50



## Endo-Flex Suction Devices

Special probes for laryngoscopy, pharyngoscopy and bronchoscopy procedures. All 12ch tubes with varying probe length. Prestrol vacuum control vent with small connector to fit Pennine connection tubing.

Product Code	Use	Length	CH/ Tip	Qty. Per Box
YS-3128/20	Laryngoscopy suction probe	20cm	12ch ø 4.00mm	10
YS-3128/26	Pharyngeal suction probe	26cm	12ch ø 4.00mm	10
YS-3128/40	Bronchoscopy suction probe	40cm	12ch ø 4.00mm	10
YS-3128/60	Bronchoscopy suction probe	60cm	12ch ø 4.00mm	10



## Sigmoid-Flex Suction Devices

An Endoflex but with standard non 'VC' vacuum control male connector. 20ch tube for Sigmoidoscopy/Oesophagoscopy.

Product Code	Use	Length	CH/ Tip	Qty. Per Box
SP-5800	Sigmoidoscopy suction probe	60cm	20ch ø 6.66mm	10




- Malleable shaft can be bent to desired position
- Longer probe
- Longer, non-slip handle

## Frazier Probes

A single-use probe manufactured from anodized aluminium for extra flexibility. With 'VC' vacuum control hole for precise control over the amount of suctioning. Used in open surgical procedures where delicate suctioning is required.

Product Code	Use	Length	CH/ Tip	Qty. Per Box
MF-1108	For neuro surgery and other intricate suction procedures	230mm	8ch	50
MF-1110		230mm	10ch	50
MF-1112		230mm	12ch	50

# 14.



## Easyflow

Filtered suction system for orthopaedic surgery.

Product Code	Use	Contents	Qty. Per Box
EFS-1111	Orthopaedic Surgery	7mm/3m Flexi-Rib Tubing - 2 curved probes - 1 straight probe - Spare filter	10

- Suitable for removal of large bone fragments during orthopaedic surgery
- Suitable for hip arthroplasty
- Eliminates blockages
- Suitable for perioperative blood collection
- Rapid filter exchange whilst maintaining suction
- Reduces risk of tissue damage and post operative infections and complications
- Ergonomic handle design ensures comfort and control
- Reversible suction tip provides extra length for intra-femoral suction


Also available with 4.5m length tubing, code EFS-1145.



## Karman Curettes

Clear, flexible curettes for uterine aspiration. Available with connector to fit to Pennine rotary swivel handle.

Product Code	Diam. mm	Qty. Per Box
AC-3004	4mm	50
AC-3005	5mm	50
AC-3006	6mm	50
AC-3007	7mm	50
AC-3008	8mm	50
AC-3009	9mm	50
AC-3010	10mm	50
AC-3011	11mm	50
AC-3012	12mm	50



## Aspiration Curettes

Aspiration Curettes are also available in boxes of 50, simply add '50' suffix to product code.


Clear, rigid plastic curettes for uterine aspiration. Moulded markings for size and correct tip orientation. Seven sizes available in either curved or straight designs. Connects to Pennine rotary swivel handle and wide bore, high vacuum tubing. Available as sets pre-attached to 2.0m thick walled aspiration tubing.

Product Code	Type	Qty. Per Box
Curved Tip		
KC-1007	7mm	10
KC-1008	8mm	10
KC-1009	9mm	10
KC-1010	10mm	10
KC-1011	11mm	10
KC-1012	12mm	10
KC-1014	14mm	10
Straight Tip		
KC-1507	7mm	10
KC-1508	8mm	10
KC-1509	9mm	10
KC-1510	10mm	10
KC-1511	11mm	10
KC-1512	12mm	10
KC-1514	14mm	10

Product Code	Type	Qty. Per Box
Curved Tip 2m Sets		
KCS-2007	7mm	10
KCS-2008	8mm	10
KCS-2009	9mm	10
KCS-2010	10mm	10
KCS-2011	11mm	10
KCS-2012	12mm	10
KCS-2014	14mm	10

Product Code	Type	Qty. Per Box
Straight Tip 2m Sets		
KCS-2507	7mm	10
KCS-2508	8mm	10
KCS-2509	9mm	10
KCS-2510	10mm	10
KCS-2511	11mm	10
KCS-2512	12mm	10
KCS-2514	14mm	10

# 16.



## Wide Bore Aspiration Tubing (Flexi-Cuff)

Product Code	Type	Qty. Per Box
KHS-3000	2m wide bore + rotary handle	10
KHS-3030	3m wide bore + rotary handle	5

Wide bore tubing without swivel handle and flexi cuff for procedures requiring high pressure suction

Product Code	Type	Qty. Per Box
AT-2000	2m wide bore	10
AT-3000	3m wide bore	5

- Sterile, double wrapped individual peel pouches
- Leak proof seals
- Easily connects to curettes and Karman catheters
- Loose male connector in each pack for connection to suction machines where required
- Tubing made of PVC
- Flexi-Cuff provides extra flexibility at a crucial point offering ease of use and improved technique



## Hy-G-Vac Suction Collection System

A complete range of compatible Pennine products designed for collection of surgical suction waste during suction procedures.

The core of the suction collection system is the Pennine Hy-G-Vac collection bottle:

- A safe, non-shatter, polypropylene bottle
- Non-removable leak-proof lid
- Quick and easy to assemble and use
- One product for all applications
- Hydrophobic block filter in the lid
- 2 litre capacity
- Safe transport after use including safe incineration properties
- HY-G-Sorb gelling agent to solidify bottle contents before disposal
- Range of trolleys, brackets and other devices to adapt and fix to existing hospital equipment

Product Code	Description	Qty. Per Box
HGV-2000	Hy-G-Vac suction collection bottle 2 litre capacity	20
HGV-2010	Link (tandem) connecting tube, F/F, 30cm, 6mm	100
HGV-2020	Specimen collection adaptor	10
HGV-2040	Ring holder for bracket or rail attachment	1
HGV-2050	HY-G-SORB gelling agent (20gm to gel approx. 1000ml water)	20
HGV-2052	HY-G-SORB gelling agent (4.0kg container)	1





## Sharps Pad

For the containment of sharps prior to disposal. Enables used sharps to be counted, controlled and disposed of safely.

Product Code	Description	Length/Size	Qty. Per Box
AP-1000	Sharps Pad	120mm x 70mm	100

- Sterile and disposable
- Adhesive sharps disposal pad
- Sealable



## Tip Cleaner

For the removal and cleaning of cauterised tissue from a diathermy pencil.

Product Code	Description	Length/Size	Qty. Per Box
TC-1000	Tip Cleaner	45mm x 45mm	300
TC-1000/100	Tip Cleaner	45mm x 45mm	100

- Sterile and disposable
- Adhesive backing for attachment to drape or trolley cover




## Diathermy Pencil

Disposable, mono-polar electro-surgical pencil for soft tissue surgery.

Product Code	Description	Qty. Per Box
EP-0001	Diathermy Pencil	20

- Standard switch button
- 3m connecting cable
- 2.5cm removable blade
- Supplied in sterile peel pouches



## Smoke Evacuation Diathermy Pencil Set

Pennine Smoke Evacuation Diathermy Pencil Sets are specifically designed to safely and effectively remove the smoke plume created by the electro-surgical pencil during surgery

Product Code	Description	Length/Size	Qty. Per Box
SEP-1000	Smoke Evacuation Diathermy Pencil	6mm ø 3m Length	10

- Pre attached to a button switch Diathermy Pencil, 3 metre cable and 3 pin plug
- Removes toxic gases, human viruses and obnoxious smells
- Ergonomic design of tubing allows for optimum manoeuvrability to prevent discomfort
- Standard connector compatible with most electro-surgical generators and smoke evacuation systems

# 18.




- Moulded PVC ring
- Supplied sterile, double wrapped
- Two strips of adhesive tape

## Camera Drape

Clear drape, to be pulled over a video camera/light cable during preparation of the operating theatre.

Product Code	Description	Diameter	Length/Size	Qty. Per Box
TP-0808/S	Camera Drape	15cm	2.4m	20



- Two bonnet covers for the x-ray tank and image intensifier head
- Integral clips for secure fit
- Brand specific for Philips and Siemens machines

## C-Arm Drape\*

Pennine 3 piece C-Arm drapes are designed to cover x-ray image intensifier equipment and protect the sterile field during surgical procedures

Product Code	C-Arm Type	Clips	Drape Dimensions (cm)	Bonnet Cover Dimensions (cm)	Qty. Per Box
LAP0011	Siemens	6 x Large Clip Fittings	230 x 40	90 x 76	24
LAP0012	Philips Libra/ Philips BV29	9 x Small Clip Fittings	205 x 31.5	74 x 78	24
LAP0013	Philips BV25	6 x Small Clip Fittings	194.5 x 31.5	74 x 78	24

\*Export only.

**Pennine**

City Gate, London Road, Derby, DE24 8WY

Tel: +44 (0)1332 794880

[marketing@penninehealthcare.co.uk](mailto:marketing@penninehealthcare.co.uk)

[www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)