

# 3M Scotchcast Plus Casting Tape



3M™ Scotchcast™ Plus Casting Tape combines the benefits of a fiberglass casting tape with the handling ease of plaster to provide orthopedic surgeons and casting personnel with the casting characteristics they need. Scotchcast Plus casting tape is available in 12 attractive, patient pleasing colors.

- **Simple application.** Quick and easy to apply without mess.
- **Water activated resin.** Gives control of set time for both basic and speciality casting applications.
- **Tack free.** Gives ease of roll unwind during application and enhanced glove movement during molding, thus improving lamination.
- **Multi-directional extensibility.** Scotchcast Plus casting tape conforms well to contours, improving cast fit and appearance.
- **Durable.** Strong, water-resistant and hard wearing.
- **Lightweight.** The fiberglass construction of Scotchcast Plus tape allows for ease of mobility.
- **Porous.** Scotchcast Plus casting tapes allow the skin to breathe.
- **Improved finish.** Scotchcast Plus tape has an improved finish to help prevent scratching skin or snagging clothes.
- **Patient pleasing colors.** Available in colors that will brighten every patient's day.

**3M** *Innovation*

# 3M Scotchcast Plus Casting Tape



## Scotchcast Plus Casting Tape (white)

| Cat No | Width          | Length       | Color | Box Qty | Case Qty |
|--------|----------------|--------------|-------|---------|----------|
| 82001  | 1in - (2.5 cm) | 2yd - (1.8m) | White | 10      | 1        |
| 82002  | 2in - (5cm)    | 4yd - (3.6m) | White | 10      | 1        |
| 82003  | 3in - (7.5cm)  | 4yd - (3.6m) | White | 10      | 1        |
| 82004  | 4in - (10cm)   | 4yd - (3.6m) | White | 10      | 1        |
| 82005  | 5in - (12.5cm) | 4yd - (3.6m) | White | 10      | 1        |

## Scotchcast Plus Casting Tape (colored)

**IMPORTANT** When ordering Scotchcast Plus colored casting tape, please ensure that you quote the correct code number. Your chosen roll width must replace the asterisk (\*) when ordering. For example, 8200\*L would become 82003L for 3in Light Blue material and 4in Red material would be 82004R.

|        |                                   |              |                  |    |   |
|--------|-----------------------------------|--------------|------------------|----|---|
| 8200*L | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Light Blue       | 10 | 1 |
| 8200*R | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Red              | 10 | 1 |
| 8200*G | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Green            | 10 | 1 |
| 8200*B | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Blue             | 10 | 1 |
| 8200*Y | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Yellow           | 10 | 1 |
| 8200*J | 2, 3, 4in - (5, 7.5, 10cm)        | 4yd - (3.6m) | Glow-in-the-Dark | 10 | 1 |
| 8200*U | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Purple           | 10 | 1 |
| 8200*E | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Grey             | 10 | 1 |
| 8200*A | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Black            | 10 | 1 |
| 8200*M | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Maroon           | 10 | 1 |
| 8200*V | 2, 3, 4in - (5, 7.5, 10cm)        | 4yd - (3.6m) | Bright Green     | 10 | 1 |
| 8200*W | 2, 3, 4in - (5, 7.5, 10cm)        | 4yd - (3.6m) | Bright Orange    | 10 | 1 |
| 8200*X | 2, 3, 4, 5in - (5, 7.5, 10, 12cm) | 4yd - (3.6m) | Bright Pink      | 10 | 1 |

## Spectra Pack - 2 rolls each black, purple, red, blue, light blue

|        |               |              |  |    |   |
|--------|---------------|--------------|--|----|---|
| 82002K | 2in - (5cm)   | 4yd - (3.6m) |  | 10 | 1 |
| 82003K | 3in - (7.5cm) | 4yd - (3.6m) |  | 10 | 1 |
| 82004K | 4in - (10cm)  | 4yd - (3.6m) |  | 10 | 1 |

## Bright Pack - 5 rolls bright pink, 3 rolls bright green, 2 rolls bright orange

|        |               |              |  |    |   |
|--------|---------------|--------------|--|----|---|
| 82002H | 2in - (5cm)   | 4yd - (3.6m) |  | 10 | 1 |
| 82003H | 3in - (7.5cm) | 4yd - (3.6m) |  | 10 | 1 |
| 82004H | 4in - (10cm)  | 4yd - (3.6m) |  | 10 | 1 |

## Application Instructions



1. Apply stockinet and cast padding.



2. Remove roll from foil pouch and dip in cool water, squeezing once. This will allow four to five minutes curing time.



3. Apply in spiral motion, over-lapping by half to two thirds on each turn.



4. Mold to aid lamination, comfort and cast strength. Weight bearing in 30 minutes.

# 3M

### Health Care

3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000  
USA  
1 800 228-3957  
healthcare@3M.com  
www.3M.com/healthcare

### 3M Canada

Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1 800 563-2921



40% Pre-consumer waste paper  
10% Post-consumer waste paper

3M and Scotchcast are trademarks of 3M.

H.I. 4509

Printed in U.S.A.

All rights reserved.

© 3M 1999 70-2008-7150-0(691)ii