

VACUETTE® BLOOD COLLECTION TUBES

ITEM N° 455071



Characterisation

VACUETTE® TUBE 8 ml CAT Serum Separator Clot Activator 16x100 red cap-yellow ring, non-ridged

Intended Use

VACUETTE® Blood Collection Tubes, Holders and Needles are used together as a system for the collection of venous blood. VACUETTE® tubes are used to collect, transport, store and process blood for testing serum, plasma or whole blood in the clinical laboratory and are for professional use.

Description

Dimension	16x100 mm
Capacity	8 ml
Cap colour	red
Ring colour	yellow

Materials

PET Polyethylene terephthalate	Tube
Paper	Label
PE Polyethylene	Cap
PP Polypropylene	Ring
Bromobutyl rubber, non-latex	Rubber stopper
Separator	Separator gel
Latex	/ none used
Phthalates	/ none used

Additive

Blood Clot Activator

Manufacturer

Legal Manufacturer:
Greiner Bio-One GmbH
Bad Haller Str. 32
A-4550 Kremsmünster

Note: Printed copies will not systematically be updated. Valid at the date of issue. Information is subject to change.

Sterilisation

Sterile interior: SAL 10^{-6} (SAL = Sterility Assurance Level)
Standard: ISO 11137

Directive

DIRECTIVE 98/79/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices Other devices / General IVD

Product Storage

Store tubes at 4-25°C (40-77°F).
NOTE: Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the tube quality (i.e. vacuum loss, drying out of liquid additives, colouring, etc.)

Shelf Life

18 months

Disposal










1. The general hygiene guidelines and legal regulations for the proper disposal of infectious material should be considered and followed.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection tubes must be disposed of in suitable biohazard disposal containers, which can then be autoclaved and incinerated afterwards.
4. Disposal should take place in an appropriate incineration facility or through autoclaving (steam sterilisation).

Packaging

24 racks of 50 pcs., 1200 pcs. in total

Label information

According to ISO 15223-1, for use in the labelling of medical devices.

Symbol	Description	Symbol	Description
	Manufacturer		Temperature limit
	Use-by date		Do not re-use
	Batch code		Consult instructions for use
	Catalogue number		In vitro diagnostic medical device
	Sterilized using irradiation		

PRODUCT SPECIFICATION

Labelling

VACUETTE 8 ml yyyy-mm-dd greinerbio-one

Surname _____
 First Name _____
 Pat.No. _____
 DOB _____ Ward _____
 Date _____ Time _____ Sig. _____

455071 Ayyymm---
 CAT Serum Sep Clot Activator

VACUETTE TUBE

8 ml CAT Serum Sep Clot Activator
 16x100 50 pcs. non-ridged

REF **455071**

LOT Ayyymm---

4.08 09:05 yyyy-mm-dd

GS1 UDI

Greiner Bio-One GmbH
 Bad Haller Str. 32,
 4550 Kremsmünster, Austria Made in Austria

VACUETTE TUBE

8 ml CAT Serum Sep Clot Activator
 16x100 red cap-yellow ring, non-ridged
 24 racks of 50 pcs., 1200 pcs. in total

REF **455071**

LOT Ayyymm---

yyyy-mm-dd

Prod. 0000000 Cart. 0000

GS1 UDI

Greiner Bio-One GmbH
 Bad Haller Str. 32,
 4550 Kremsmünster,
 Austria Made in Austria

Note: Printed copies will not systematically be updated. Valid at the date of issue. Information is subject to change.