



| | |
|-----------------------------------|---|
| Description: | Absorbent dressing pad |
| Final dimensions: | 10 cm x10 cm |
| Composition: | Cellulose, Cellulose pulp, Polypropylene nonwoven, SSMMS blue PP NW (5 layers). |
| Construction: | Absorb. pads/Pensos |
| Latex: | Free |
| PVC: | Free |
| Re use: | No |
| Width variation: | ± 5 mm |
| Length variation: | ±10 mm |
| Colour: | white/ blue |
| Sterilized by: | Ethylene Oxide |
| Shelf life: | 60 months |
| CFU/g | N.A. |
| Free swell absor. capacity | >65 g/100cm ² |
| Absorption Cap. (ISO 1194) | Approx 12,5 ml |
| Absorption EN 1644-2 | Approx. 8 g/g |



| | |
|------------------------------|--|
| Nr and Revision Date: | 5 -- 1996/10/24 |
| CE Class: | I Steril |
| CE mark: | CE0120 |
| Directive: | MDD 93/42/ EEC amended by 2007/47/EC |
| Applicable norms: | EN ISO 13485:2016 ; EN ISO 11607-1:2009 ; EN ISO 11607-2:2006 ; EN ISO 11737-1:2006 ; EN ISO 11737-2:2009 ; EN 556-1:2001 ; EN ISO 13726-2:2002 ; EN ISO 14971:2012 ; EN ISO 10993-1/AC:2010 ; EN ISO 10993-7:2009 ; EN 1422:2014 ; EN ISO 11737-1:2006/AC:2009 ; EN 556-1:2001/AC:2006 ; EN ISO 11135:2014 ; ISO 16142-1:2016 ; EN ISO 15 223-1:2016 ; |

| Level | Type | Dimension | Stamp | Labels |
|-------|---------------------------|---------------|---------|---------|
| 1 | 1 Flexible peel-open pack | 134x140mm | stamped | |
| 2 | 60 Folding box | 142x315x150mm | | labeled |
| 3 | 600 Carton box | 780x285x315mm | | labeled |

Device can be subject to cumulative single, multiple or repeated use or contact up to 30 days (prolonged; type B) with intact, breached or compromised body surface.

| | |
|-----------------------------|--|
| Storage instructions | Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply |
| Waste instructions | No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste. |
| Intended use | Primary dressing intended for use on moderately to heavily exuding wounds and open drains. Universal use as absorbent pad in general patient care. Ideal for situations where a great capacity of absorption is needed, such as bleedings and chronicle wounds with moderate to high exudation. Secondary dressings for high exuding wounds. |

| | | |
|--|---|-----|
| WASTE AND ENVIRONMENT REGULATIONS | Weight (gr) of packing material used in 100 pcs | 544 |
|--|---|-----|

ISSUED BY: Quality Management Department, 2019/08/08