




Description:	Cellulose wadding
Final dimensions:	40 cm x60 cm
Composition:	Cellulose.
Latex:	Free
PVC:	Free
Re use:	No
Sterilizable by:	Ethylene oxide
Colour:	white
Shelf life:	120 months
CFU/g	<10 ²
Size variation	±10 %



Nr and Revision Date:	6 -- 2019/01/08
Directive:	Not applicable
Applicable norms:	EN ISO 13485:2016 ; EN ISO 11737-1:2006 ; ENISO11737-1:2006/AC:2009 ;

Important Notes:	5 Kg per plastic bag
-------------------------	----------------------

Level	Type	Dimension	Stamp	Labels
1	5 Plastic bag			labeled
2	20 Carton box	590 x 390 x 408		labeled

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has good stability and therefore, normal conditions of medical devices storage shall apply.
Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
Intended use	To cover surfaces

WASTE AND ENVIRONMENT REGULATIONS	Weight (gr) of packing material used in 100 pcs	4175
--	---	------

ISSUED BY: Quality Management Department, 2019/08/08

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.