

Levins Tube Products

Product Data Sheet

Description	Levins Tube, Graduated, X-Ray with open tip, four side eyes and colour coded funnel connector.
Intended Use	Levins tube is used for gastro-intestinal feeding and aspiration of intestinal secretions. Levins tube is a tube that is inserted through the nose into the upper alimentary canal and is used to facilitate intestinal decompression. *Please note Pennine does not recommend non ENFit tubes for feeding in markets that have adopted the ENFit System.
Duration of Use	Short Term Use: for more than 60 minutes and less than 30 days.
Materials	Tube & Funnel Connector - PVC (DEHP Free); Radio-Opaque Line (X-Ray). Not made with DEHP Plasticizer. Not made with Natural Rubber Latex.
Product Information	Smooth outer surface of the tube to reduce chances of tissue irritation. Soft, frosted and kink resistance PVC tubing. 125cm tube with radio-opaque line provided to its full length and marked at 50, 60, and 70 cm from the tip for accurate placement of tube. Open end (terminal eye) to ensure maximum dispersion. Four lateral eyes to reduce risk of blockage. Color coded funnel end connector for easy identification of size. Tube shore hardness: 75~78A. Available in different sizes from FG 6 to FG 20. Supplied sterile, in individual peel pouches.
Quality System	EN ISO 13485:2016
EU Classification	Class IIa
Sterilised by	Ethylene Oxide
Shelf Life	5 years
Storage & Safety	Toxicity: Manufactured from biocompatible materials. Storage: Store in between 5°C to 35°C, avoid excessive heat, protect from direct sunlight and moisture.

Product Code	Tube Length	CH/Dia	Connector	Side Eyes	Tip	Box Qty
AEXLT2206	125cm	6ch ø 2.0mm	Funnel	4	Open Tip	50
AEXLT2208	125cm	8ch ø 2.7mm	Funnel	4	Open Tip	50
AEXLT2210	125cm	10ch ø 3.3mm	Funnel	4	Open Tip	50
AEXLT2212	125cm	12ch ø 4.0mm	Funnel	4	Open Tip	50
AEXLT2214	125cm	14ch ø 4.7mm	Funnel	4	Open Tip	50
AEXLT2216	125cm	16ch ø 5.3mm	Funnel	4	Open Tip	50
AEXLT2218	125cm	18ch ø 6.0mm	Funnel	4	Open Tip	40
AEXLT2220	125cm	20ch ø 6.7mm	Funnel	4	Open Tip	40

***Please note Pennine does not recommend non ENFit tubes for feeding in markets that have adopted the ENFit System.**

© 2023 Pennine Healthcare or one of its subsidiaries.
All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



Poly Medicure Limited, Plot No. 104-105, Sector-59, HSIIDC Industrial Area, Ballabhgarh, Faridabad, Haryana 121004, INDIA




Obelis S.A. Boulevard General Wahis 53,
B-1030, Brussels, Belgium. mail@obelis.net



Levins Tube Products

Product Data Sheet



*Please note Pennine does not recommend non ENFit tubes for feeding in markets that have adopted the ENFit System.

© 2023 Pennine Healthcare or one of its subsidiaries.
All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



Poly Medicure Limited, Plot No. 104-105, Sector-59, HSIIDC Industrial Area, Ballabhgarh, Faridabad, Haryana 121004, INDIA



Obelis S.A. Boulevard General Wahis 53,
B-1030, Brussels, Belgium. mail@obelis.net

