


Description:	Kocher forceps
Final dimensions:	13 cm
Composition:	Reinforced Polypropylene.
Latex:	Free
PVC:	Free
Re use:	No
Sterilizable by:	Ethylene Oxide
Colour:	blue
Shelf life:	120 months
CFU/g	<10 ²



Nr and Revision Date:	27 -- 2019/08/09
CE Class:	I
CE mark:	CE
Directive:	MDD 93/42/EEC amended 2007/47/EC
Applicable norms:	EN ISO 13485:2016 ; EN ISO 11737-1:2006 ; EN ISO 14971:2012 ; EN ISO 10993-1/AC:2010 ; ISO 1133-1:2011 ; ENISO11737-1:2006/AC:2009 ; ISO 16142-1:2016 ; EN ISO 15 223-1:2016 ;

Level	Type	Dimension	Stamp	Labels
1	50 Plastic bag			labeled
2	2000 Carton box	590x390x408mm		labeled

Important Notes:	Nr pcs/pack subjected to ± 3% variation
-------------------------	---

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of storage for medical devices shall
-----------------------------	--

Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.
---------------------------	---

Intended use	Intended to grasp and hold dressing material and other items during medical procedures.
---------------------	---

WASTE AND ENVIRONMENT REGULATIONS	Weight (gr) of packing material used in 100 pcs	65
--	---	----

ISSUED BY: Quality Management Department, 2019/08/09