

iHealth Fit

Smart body analysis scale HS2S

Compatible with iPhone and Android smartphones



Weight and body composition monitoring on iPhone and Android smartphones

Requires iOS v9.0 and + or Android version 5.0 and +

iHealth Fit features

- Light and compact body analysis scale
- Gives weight, BMI (Body Mass Index), lean, bone, muscle, water mass, body fat, Daily caloric intake and Visceral fat rating
- Luminous display (LED)
- 8 users
- Tracks data, graphs and history on the iHealth MyVitals App
- Store up to 500 measurements per user before synchronization.
- Share your data with your healthcare professional or family
- Bluetooth Sync, easy to install
- Secure and free data storage



Technical specifications:

- Connectivity: Bluetooth 4.0
- Dimensions: 306mm × 306mm × 26.6mm
- Weight: 1850g
- Battery: 4 AAA batteries 1.5V
- Range:
 - Weight: 5kg-180kg (11 lb-400 lb)
 - Body fat: 5%-65%
 - Resolution: 0,1kg / 0.1lb
- Accuracy:
 - Weight: ±1.1lb / 0.5 kg (11-110lb / 5kg-50kg); ±1% (110-400lb / 50kg-180kg)
 - Body fat: ±1%
- Operating Environment:
 - Temp: 5°C ~ 40°C, humidity: 20-85% RH
- Conditions for transport and storage:
 - Temp: -20°C ~ 55°C, humidity ≤90% RH



In the box:

- 1 iHealth Fit body analysis scale
- 1 quick start guide
- 4 AAA batteries



Track your vitals information with the free iHealth MyVitals App



Compatible with all the iHealth products:

Activity tracker, blood pressure and SpO2 monitors



www.iHealthLabs.eu

iHealthTM
Connected to my health