

Silicone Respiration Mask

REF: 51.500.00.000 - 005

Unit: 1 pcs

Packing Unit: (see below)

Description

Silicone masks are used for respiration of patients. The masks can be reused and prepared.

REF	size	Size	Packi ng unit	holder
51.500.00.000	Mask for neonatal	0	50	
51.500.00.001	infant size 1, round	1	50	
51.500.00.002	infant size 2, round	2	50	
51.500.00.003	pediatric size 3	3	50	x
51.500.00.004	adult size 4	4	50	x
51.500.00.005	adult size 5	5	100	x

Material Composition

Product: mask: Silicone
holder: PC

Packing: individual packing: PE
outer case: cardboard

Sterilisation

Non-sterile