

SPECIFICATIONS

ROMED LEMON GLYCERIN SWABSTICKS

REF: LSWAB-250



Medical device conform MDR (EU) 2017/745

Soft cotton tips for gentle oral lubrication and hydration.
Pleasant tasting lemon solution.
Soothes and refreshes dry mouths of patients with restricted oral intake.

Materials

Solution:

Glycerin (60%)
Citric acid (3%)
Sodium benzoate (0.1%)
Lemon extract (< 0.1%)
Water (37% ± 1%)

Amount of impregnation solution per pack:

5 ml

Material plastic applicator:

Polypropylene

Total length:

100 mm (± 5%)

Length cotton head:

30 mm (± 5%)

Diameter cotton head:

10 mm (± 2 mm)

Cotton standard:

British Pharmacopeia 5.0

Packaging material:

Paper, Aluminum, Polyethylene

Free from pesticides and other contaminating agents. No ingredients of animal origin. Latex free.
Contains no sugar and can be used for diabetic patients.

Validity of 5 years if package is undamaged.

Packing

Three pieces in a lightproof pack, 25 packs in an inner box, 10 inner boxes in an export carton