

# Product specifications

## Arseus

Made by Ab Kronoby Fatex Oy

Quality requirements	Slabbetje 3L ronde hals, beige 38x70		
Testing parameter	Dim.	Coc	2113270
<b>Appearance, property</b>			
- Color			Beige
- Dirt			Free from dirt, lint and stains
- Neckline			Band fastening
- Pocket			Turnable pocket
<b>Components</b>			
- Tissue (TCF-pulp)			Towel Tissue 18gsm/800mm
- Tissue (TCF-pulp)			Towel Tissue 18gsm/600mm
- PE-film			PE-film 900x0,014mm Beige 2-S
- Adhesive			-
- Extras			Folded PE-film by the neckline to produce strengthening
<b>Length</b>	mm	Yes	700 ± 20mm
<b>Width</b>	mm	Yes	380 ± 10mm
<b>Length (folded)</b>	mm		195
<b>Width (folded)</b>	mm		380
<b>Weight</b>	g	Yes	14,09g ± 15%
<b>Water penetration</b>		Yes	0
<b>Bioburden</b>	CFU/g	Yes	≤1000CFU/g
<b>Packaging</b>			
PE-bag, Qty/bag	Pcs		100
PE-bag, Qty/Carton	Pcs		6
Cartons/Pallet	Pcs		48
Carton Size (Inner)	mm		387x387x280
PE-bag Size (Length x Width)	mm		300x500
<b>Pathogenic germs</b>	CFU/g	Yes	Free from pathogenic germs like staph. Aureus, Enterobacteriae, E. coli, Pseudomonas aeruginosa
<b>Check of certificate</b>		Yes	

Quality requirements			Washandje Airlaid met PE
Testing parameter	Dim.	Coc	3118000
<b>Appearance</b>			
- Color			Blue
- Dirt			Free from dirt, lint and stains
<b>Components</b>			
- ALT (TCF-pulp)			Airlaid 50g/m <sup>2</sup> 450mm
- PE-film			PE-film 440x0,014mm Blue
- Adhesive			0,204g/pcs
<b>Length</b>	mm	Yes	230 ± 10mm
<b>Width</b>	mm	Yes	155 ± 5mm
<b>Weight</b>	g	Yes	4,75g ± 15%
<b>Water penetration</b>		Yes	0
<b>Bioburden</b>	CFU/g	Yes	≤1000 CFU/g
<b>Packaging</b>			
PE-bag, Qty/bag	Pcs		50
PE-bag, Qty/Carton	Pcs		28
Cartons/Pallet	Pcs		42
Carton size (Inner)	mm		390x390x325
PE-bag Size (Length x Width)	mm		215x350
<b>Pathogenic germs</b>		Yes	Free from pathogenic germs like staph. Aureus, Enterobacteriae, E. coli, Pseudomonas aeruginosa
<b>Check of certificate</b>		Yes	