



# iHealth Air

## Connected Pulse Oxymeter P03m

Compatible with iPhone, iPad and Android smartphones



**Monitor oxygen saturation, perfusion index and pulse on iPhone, iPad and Android smartphones**

For iOS and Android version 4.0 +.

See the list of compatible devices on

[ihealthlabs.eu/support](http://ihealthlabs.eu/support)

### iHealth Air features

- Measures and instantly displays on the device and your smartphone:
  - SpO2 : Concentration of oxygen in the blood
  - Pulse : Number of heartbeats per minute
  - PI : Perfusion Index displayed on iHealth MyVitals App
- Tracks data, graphs and history on the iHealth MyVitals App
- Bluetooth 4.0 wireless synchronization
- Share your medical data with you healthcare professional or family
- Charges by micro USB
- Secure and free data storage



### Technical specifications

- Connectivity: Bluetooth 4.0 BLE
- Display: LED
- Size: 62 x 33 x 28mm
- Weight: 42g
- Power supply: battery, 3.7V Li-ion, 300mAh
- Measuring range SpO2: 70-99%
- Measuring range Pulse: 30-250bpm
- Accuracy SpO2: 80-99%, ± 2%; 70-79%, ± 3%
- Accuracy of pulse rate: ± 2bpm
- Automatic shutdown after 8 seconds
- Operating environment:  
Temp: 5 - 40 °C; Humidity <80% RH
- Storage environment:  
Temp: -20 - 55 °C; Humidity <95% RH



### In the box:

- 1 iHealth Air Pulse Oxymeter
- 1 neck strap
- 1 USB cable
- 1 quick user guide
- 1 user manual

Track your vital information with the free iHealth MyVitals App



Compatible with many other great iHealth products:  
Activity trackers, blood pressure and scales



[www.iHealthLabs.eu](http://www.iHealthLabs.eu)

**iHealth**<sup>TM</sup>  
Connected to my health