



Micropump™
MP-101+
MP Daily+

Freedom without clinical compromise

 micrel

Syringe pumps designed for portability, freedom and simplicity

Micrel MP-101+ and MPDaily+ syringe pumps offer accurate delivery, are easy to maintain with a trouble-free set-up



Simple set up



A variety of syringe sizes and brands can be used



Long battery life



Long service interval

- ✓ A clear informative display during programming and easy alarm identification
- ✓ Incorporating a dual micro-processor for extra confidence during drug delivery applications
- ✓ Suitable for hospital and home care use for pain management and other infusions

- Use any manufacturer's standard syringes (i.e BD, Terumo, BBraun etc) from 10 to 50 ml
- Event history log to store and review pump activities when needed
- Safe and convenient syringe attachment
- Can be comfortably worn with belt or shoulder holster
- Backlit LCD display for clear information on flow rate and identification of alarms

SPECIFICATIONS

MP-101+



MPDaily+



Dimensions:

170 x 61 x 32 mm

Rate set up:

mm/hr

mm/24hr

Flow rate:

1 to 99 mm/hr (1 mm increments)

1 to 99 mm/24hr (1 mm increments)

Accuracy:

±2%

Weight:

220 g including batteries (Alkaline LR03 AAAx6)

Battery life:

approx. 2 months with 50mm/day infusion

CONFIGURABLE

- Accurate bolus delivery with lock out feature
- Configurable occlusion level
- Lockable flow rate
- Simple priming procedure

The pumps have audio and visual alert modes: alarms and warnings (i.e low battery/depleted battery, end of infusion, multi point syringe dislodged detection, occlusion etc.)

ACCESSORIES

Each device includes a carrying bag, operating instructions and a set of batteries.

Soft locker box is available for 20ml & 50ml syringes.
Compact pole clamp. Single use bag for MP+ pumps.



Discover a patient-centred solution

To find out more about MP-101+ & MPDaily+ syringe pumps go to micrelmed.com