



OPERATING ROOM
Skin adhesive



derma+flex



Value Life

derma+flex

The **strength** of Butyl & the **flexibility** of Octyl

High viscosity skin adhesive

Indications

derma+flex is a topical skin adhesive for wound edges from surgical incisions and also incisions from minimally invasive surgery, and simple, thoroughly cleansed, trauma-induced lacerations.

Benefits for

Medical Staff

Precise wound closure thanks to this unique applicator

- Easy to apply
- Creates a microbial barrier
- Suture strength
- Highly flexible
- Suture & staple compatible
- Fast wound closure (60-90sec)
- Duo tip allows for both precise and broad application
- No sharps risks

Patient

Increased patient comfort

- Waterproof, shower friendly formula
- Painless application
- Does not require a return visit for removal
- Minimal exothermic effect (it won't burn or sting)
- Good aesthetic results

Hospital


Logistic & economical benefits

- Lower cost
- Room temperature storage
- Safe and reliable
- Shelf life: 24 months

Features

Dual Use applicator which is already **pre-assembled** on the tube.

Unique formula of
60% 2-octyl cyanoacrylate: strength
40% n-butyl cyanoacrylate: flexibility



Product	Code
derma+flex 0,7g	VCHLB60406V

CARDIOVASCULAR SURGERY

For further information, please contact: marketingbenelux@vygon.com

The specifications given in this brochure are for information only
and are not, under any circumstances, of a contractual nature.

SA Vygon NV • Haachtsesteenweg 1650 Chaussée de Haecht • 1130 BRUSSELS • BELGIUM
Tel: +32 (0)2 706 09 50 • www.vygon.be
Vygon Nederland B.V. • Kerkhofstraat 21 • 5554 HG VALKENSWAARD • NEDERLAND
Tel: +31 (0)40 208 93 80 • www.vygon.nl
VY21010 •  Vygon Benelux

