

1541510 - OP-C2L colposcope avec LED intégré



Déscription

The OP-C2L colposcope combines the standard qualities of the OPTOMIC colposcopes with the latest LED technology illumination. The LED, integrated in the colposcope head, offers white and uniform lighting, perfect for carrying out the daily work of the specialist.

- LED illumination standard
- 3 different magnifications
- Light and compact
- Inclined binocular head
- Optional video system

Spécifications

Marque	Optomic
sku	1541510
Réglementation Européenne	MDR - 2017/745/EU - Classe I
Type d'emballage	Pièce