

## 1541510 - OP-C2L colposcope avec LED intégré



### Déscription

The OP-C2L colposcope combines the standard qualities of the OPTOMIC colposcopes with the latest LED technology illumination. The LED, integrated in the colposcope head, offers white and uniform lighting, perfect for carrying out the daily work of the specialist.

- LED illumination standard
- 3 different magnifications
- Light and compact
- Inclined binocular head
- Optional video system

### Spécifications

<b>Marque</b>	Optomic
<b>sku</b>	1541510
<b>Réglementation Européenne</b>	MDR - 2017/745/EU - Classe I
<b>Type d'emballage</b>	Pièce