

1013437 - Eucerin pH5 lotion nettoyante - pompe - 1000 ml - 1 pc



## Déscription

La lotion lavante Eucerin pH5 est un produit de lavage quotidien du visage et du corps qui ne contient pas de savon alcalin. Les tensioactifs extra-doux nettoient en douceur mais efficacement les peaux sèches et sensibles. Le système unique Eucerin pH Balance donne au nettoyant la valeur pH optimale pour protéger les défenses naturelles de la peau, maintenir son élasticité et empêcher son dessèchement. Il convient également aux bébés et aux enfants.

## Spécifications

<b>Marque</b>	Beiersdorf
<b>Tube</b>	Non
<b>Système de pompage</b>	Oui
<b>Ingrédients</b>	AQUA, SODIUM COCOAMPHOACETATE, SODIUM MYRETH SULFATE, LAURYL GLUCOSIDE, CITRIC ACID, PEG-40 HYDROGENATED CASTOR OIL, PEG-200 HYDROGENATED GLYCERYL PALMATE, COCO GLUCOSIDE, GLYCOL DISTEARATE, GLYCERINE, POLYQUATERNIUM-10, DIAMMONIUM CITRATE, SODIUM BENZOATE, SODIUM SALICYLATE, LINALOOL, BUTYLPHENYL METHYLPROPIONAL, ALPHA-ISOMETHYL IONONE, HEXYL CINNAMAL, PARFUM.
<b>sku</b>	1013437
<b>Contenu</b>	1000 ml
<b>Forme physique</b>	Savon liquide
<b>Type d'emballage</b>	Pièce