

1013437 - Eucerin pH5 lotion nettoyante - pompe - 1000 ml - 1 pc



Déscription

La lotion lavante Eucerin pH5 est un produit de lavage quotidien du visage et du corps qui ne contient pas de savon alcalin. Les tensioactifs extra-doux nettoient en douceur mais efficacement les peaux sèches et sensibles. Le système unique Eucerin pH Balance donne au nettoyant la valeur pH optimale pour protéger les défenses naturelles de la peau, maintenir son élasticité et empêcher son dessèchement. Il convient également aux bébés et aux enfants.

Spécifications

Marque	Beiersdorf
Tube	Non
Système de pompage	Oui
Ingrédients	AQUA, SODIUM COCOAMPHOACETATE, SODIUM MYRETH SULFATE, LAURYL GLUCOSIDE, CITRIC ACID, PEG-40 HYDROGENATED CASTOR OIL, PEG-200 HYDROGENATED GLYCERYL PALMATE, COCO GLUCOSIDE, GLYCOL DISTEARATE, GLYCERINE, POLYQUATERNIUM-10, DIAMMONIUM CITRATE, SODIUM BENZOATE, SODIUM SALICYLATE, LINALOOL, BUTYLPHENYL METHYLPROPIONAL, ALPHA-ISOMETHYL IONONE, HEXYL CINNAMAL, PARFUM.
sku	1013437
Contenu	1000 ml
Forme physique	Savon liquide
Type d'emballage	Pièce