

## 1501401 - LTV 1200 ventilateur



### Déscription

Transport ventilator Simple to use, advanced ventilation in a small lightweight package, that's the LTV 1200 ventilator. Use the LTV 1200 in the ICU, emergency department or for patient transport with confidence that you are delivering the highest quality in patient ventilation care and comfort. Set up initial ventilation quickly on the LTV 1200 ventilator with presets for Infant, Pediatric or Adult patients using either invasive or noninvasive modes. Refine the settings easily using the touchturntouch interface with the LED illuminated front panel. For patients 5 kg and above, the LTV 1200 ventilator provides a wide range of ventilation therapies to meet demanding patient needs, including Volume Control, Pressure Control, Pressure Support and Spontaneous breath types. The ultra sensitive flow trigger combined with the Spontaneous Breathing Trial function makes weaning easier and more efficient for both the patient and clinician. The LTV 1200 ventilator gives you the freedom of conveniently portable ICU ventilation.

### Spécifications

<b>Marque</b>	Vyaire
<b>sku</b>	1501401
<b>Réglementation Européenne</b>	MDD - 93/42/EEC - Classe IIa
<b>Type d'emballage</b>	Pièce