

1603833 - B. Braun kit d'anesthésie épidurale - 1 pc



Spécifications

Marque	Braun
sku	1603833
Contenu des sets stériles	Perifix® LS LOR seringue 8 ml: 1pc, Omnifix® seringue LS 10 ml: 1pc, Omnifix® LS seringue 5 ml: 1pc, Omnifix® LS seringue 3 ml: 1pc, Sterican® canule 26G, 0,45 x 25 mm, brun: 1pc, Sterican® canule 21G, 0,8 x 40 mm, vert: 1pc, Sterican® canule G18, 1,2 x 40 mm, rose: 1pc, compresses de gauze 7,5 x 7,5 cm, 8-plis: 5pcs, swabsticks: 2pcs, drap patients autocollant, 75 x 90 cm, bleu: 1pc, Plateau 20 x 14,8 x 4 cm, vert: 1pc, plateau rectangulaire, 14,8 x 7,5 cm, 265 ml, blanc: 1pc, champs vert, 75 x 90 cm: 1pc
	Non
Réglementation Européenne	MDR - 2017/745/EU - Classe III
Type d'emballage	Set