

## 1550123 - Schiller Trolley X1 avec fixation pour Cardiovit FT-1



### Déscription

The Trolley X1 is part of a new family of SCHILLER trolleys. With its modern and ergonomic design, the X1 guarantees maximum mobility and steering control, as well as functional cable management through its aluminum column. Equipped with a specific velcro-free adhesive mounting bracket for the FT-1.

### Spécifications

<b>Marque</b>	Schiller
<b>sku</b>	1550123
<b>EKG catégorie</b>	Divers ECG
<b>Type d'emballage</b>	Pièce