

## 1021249 - Eucerin pH5 lotion nettoyante sans savon avec pompe - 400 ml - 1 pc



### Déscription

La lotion lavante Eucerin pH5 est un produit de lavage quotidien du visage et du corps qui ne contient pas de savon alcalin. Les tensioactifs extra-doux nettoient en douceur mais efficacement les peaux sèches et sensibles. Le système unique Eucerin pH Balance donne au nettoyant la valeur pH optimale pour protéger les défenses naturelles de la peau, maintenir son élasticité et empêcher son dessèchement. Il convient également aux bébés et aux enfants.

## Spécifications

<b>Marque</b>	Beiersdorf
<b>Ingrédients</b>	Aqua, sodium cocoamphoacetate, sodium myreth sulfate, lauryl glucoside, sodium chloride, citric acid, PEG-40 hydrogenated castor oil, PEG-200 hydrogenated glyceryl palmate, coco glucoside, glycol distearate, glycérine, polyquaternium-10, diammonium citrate, sodium benzoate, sodium salicylate, benzyl alcohol, linalool, butylphenyl methylpropional, alpha-isomethyl ionone, hexyl cinnamal, benzyl salicylate, parfum
<b>sku</b>	1021249
<b>Contenu</b>	400 ml
<b>Forme physique</b>	Savon liquide
<b>Odeur</b>	Parfumée
<b>Type d'emballage</b>	Pièce