

1541510 - OP-C2L colposcoop met geïntegreerde led lichtbron



Beschrijving

The OP-C2L colposcope combines the standard qualities of the OPTOMIC colposcopes with the latest LED technology illumination. The LED, integrated in the colposcope head, offers white and uniform lighting, perfect for carrying out the daily work of the specialist.

- LED illumination standard
- 3 different magnifications
- Light and compact
- Inclined binocular head
- Optional video system

Specificaties

Merk	Optomic
sku	1541510
Europese Regelgeving	MDR - 2017/745/EU - Klasse I
Type verpakking	Stuk