

**1543939 - Renol pu sloop armkussen - h 15 cm x b 20 cm x l 35 cm
- 1 st**



Beschrijving

Afneembare overtrek in katoen en wasbaar in de machine.

Specificaties

Merk	Renol
SKU leverancier	RCM16099
Normeringen	EN ISO 13485:2016; 14971:2012; 10993:2019; 10993-5:2009; 1041:2008; 15223-1:2016; 62366:2008
Samenstelling	PU
sku	1543939
Kleur sloop	Blauw
Materiaal	PU
Merk positioneringsmateriaal	Renol
Wasbaar	60 °C
Europese Regelgeving	MDR - 2017/745/EU - Klasse I
Type verpakking	Stuk