

1543930 - Renol positioneringskussen lang - katoen - diameter 26 x 35 cm x l 280 cm - 1 st



Specificaties

Merk	Renol
Product afmetingen	B 35 cm, L 280 cm
SKU leverancier	RCM12233
Normeringen	EN ISO 13485:2016; 14971:2012; 10993:2019; 10993-5:2009; 1041:2008; 15223-1:2016; 62366:2008
Samenstelling	Katoen
sku	1543930
Kleur sloop	Blauw
Materiaal	Katoen
Merk positioneringsmateriaal	Renol
Wasbaar	nvt
Europese Regelgeving	MDD - 93/42/EEC - Klasse I
Type verpakking	Stuk