

## 1541858 - Ergoselect 50P



### Beschrijving

The cost-effective ergometer for exercise ECG tests. A low step-through frame, freely definable exercise protocols and universal interfaces for connecting to all commercially available PC-based ECG systems and ECG recorders set this ergometer apart.

Accessories: Horizontal saddle adjustment, Racing saddle, Quick-release adapter for standard saddle, Adjustable pedal cranks.

### Specificaties

<b>Merk</b>	Ergoline
<b>sku</b>	1541858
<b>Type verpakking</b>	Stuk